

Laser induced deflection (LID) – absolute photo-thermal absorption measurement in optical coatings and materials

Dr. Christian Mühlig



Leibniz Institute of Photonic Technology, Jena, Germany

Outline

- Principle of **L**aser **I**nduced **D**eflection (LID) technique
- Calibration procedure
- Setup
- Measurement concepts
- Experimental results
- Summary

Laser induced deflection (LID)




Photo-thermal technique with *transversal* probe beam guiding



- Probe beam deflection is a *direct measure* for the absorbed mean laser power
- **Main challenge:** How to calibrate for *reliable* absolute absorption data?


Calibration (electrical)

- Absorption induced heat is simulated by electrical heating (adaption of calorimetry's calibration to a photo-thermal technique)
- **Absolute data without simulations and knowledge of material's photo-thermal parameters**
- Unique feature for photo-thermal techniques



- Linear behavior over multiple orders of magnitude of absorbed power
- Energy balance measurements ($T+R+A+S=1$) approve that no systematic error exists




Setup




Example for sensitivity (fused silica): $10\mu\text{W}$ absorbed power

Measurement concepts


Different concepts for different requirements



Measurement concepts





„Sandwich“ concept





Separating materials for pump and probe beam guiding allows for:

- ✓ Measuring small samples like NLO and laser crystals
- ✓ Measuring samples without side-face polishing
- ✓ Strong sensitivity increase for materials with low photo-thermal response by choosing appropriate sandwich tile materials
- ✓ Measuring samples which are not transparent for VIS probe beam wavelengths
- ✓ Simplified electrical calibration for bulk absorption




Measurement concepts



Concepts for particular „coating“ substrates



- (Large) rectangular samples (e.g. 20x20x6mm³)



- $\varnothing=1/2''$, $d=1/4''$
- (Small) rectangular samples





- $\varnothing=1...2''$, $d=1\text{mm}$




Experimental results


Measurement example



Experimental results




Measurement example




- Irradiation induced absorption increase to steady state value
- MgF_2 part of the AR seems to be responsible

Experimental results




Measurement example







- Irradiation induced absorption increase to steady state value
- MgF_2 part of the AR seems to be responsible


Experimental results



Absorption as function of particular impurities or deposition parameters





Experimental results




Thin film thickness series to separate coating and interface absorption

- Fluoridic single layers (optical thickness: multiple of $\lambda/4$) on thin CaF_2 disc, attached to fused silica cube for sensitivity enhancement
- Total absorption vs. film thickness allows to separate between interface and layer absorption
- Requirements: identical CaF_2 discs, interface and layer properties



Experimental results





Effect of IBF treatment on CaF₂ surface absorption at 193nm

Sample	Treatment	Absorption (% per surface @193nm)
# 1	Pitch polishing (standard)	0.08 ± 0.008
# 2	Pitch polishing (standard)	0.083 ± 0.01
# 5	IBF (standard parameters)	0.3 ± 0.035
# 6	IBF (double neutralization)	0.39 ± 0.03
# 7	IBF (double treatment time)	0.445 ± 0.035
# 8	IBF (step-like treatment)	0.355 ± 0.03
# 9	IBF (optimized removal)	0.365 ± 0.035
# 10	IBF (preserve treatment)	0.355 ± 0.03

IBF treatment



Enhanced surface absorption




Sample	Absorption (% per surface)	Absorption after coating process (%)		Absorption decrease after coating, taking into account an AR absorption of ~0.0003 (%)	
		AR coated side	Uncoated rear side	AR coated side	Uncoated rear side
# 1	0.08 ± 0.008	0.121 ± 0.016			
# 2	0.083 ± 0.01	0.125 ± 0.015			
# 5	0.3 ± 0.035	0.12 ± 0.012	0.23 ± 0.025	0.21	0.07
# 7	0.445 ± 0.035	0.118 ± 0.01	0.225 ± 0.025	0.357	0.22
# 9	0.365 ± 0.035	0.125 ± 0.012	0.25 ± 0.025	0.27	0.115
# 10	0.355 ± 0.03	0.135 ± 0.02	0.27 ± 0.025	0.25	0.085



After one-side AR coating (LaF₃/MgF₂):

- Total absorption of AR coated sides are virtually identical for IBF and pitch polished samples.
- Uncoated rear side shows higher absorption than AR coated side.
- Proposed annealing mechanisms:
 - High deposition temperature
 - Fluorine supply by coating materials

Experimental results



Combining absorption and CRD measurements



LID-measurement:
 $A(532\text{nm}) = 18.8 \text{ ppm}$

From multiple CRD-measurements:
 $(T+A+S) = 113 \text{ ppm}$

Summary

- **Laser Induced Deflection (LID)** allows for highly sensitive direct absorption measurements by using the photo-thermal effect
- Unique electrical calibration procedure provides verified absolute absorption data without knowledge of material parameters/simulations
- Table-top setup with access to all wavelengths if required (laser) power is available (no restrictions to laser beam quality)
- Variety of measurement concepts for different requirements
- Combination with additional measurement techniques (CRD, fluorescence...)

Acknowledgement

- The Federal Ministry of Economics and Technology for financial support within the framework of the Confederation of Industrial Research Associations (projects: “LIDAG-OK” and “Kompakter, gepulster 266nm-Festkörperlaser hoher Leistung und Lebensdauer für die Materialbearbeitung“).
- Numerous partners for their trust in our measurement technique and ongoing new challenges.

There is also a poster on the topic for
more information
and time for discussion!